trophon[®]2

User Manual Addendum

for trophon[®]2 Plus

Software Upgrade

ENGLISH	3
DEUTSCH.....	14
FRANÇAIS	25
NEDERLANDS	36
ITALIANO	47
NORSK.....	58
SUOMI	69
DANSK	80
SVENSKA.....	91
ESPAÑOL.....	102
TÜRKÇE.....	113

trophon[®]2

User Manual Addendum

for trophon[®]2 Plus

Software Upgrade

Read this addendum together with trophon®2 user manual before operating device to determine the correct procedures.

Table of Content in this addendum lists updates to the trophon®2 User Manual.

For further information, contact your customer service representative or visit the Nanosonics website.

© 2025 Nanosonics Limited. All rights reserved.

www.nanosonics.com

Table of Contents

Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS 6
 A1.1 Labels and Symbols.....6
 A2.1 Indications for Use8

Part B – SETUP 9
 B2.6 Basic Settings Network9


Part C – CYBERSECURITY 10



GLOSSARY..... 13

Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS

The following section of the labels and symbols table has been updated to add a new feature relating to DICOM in Touch Screen Symbols.

A1.1 Labels and Symbols

	Caution		Warning
	Consult Instructions For Use		Corrosive
	Environmental Conditions: trophon Storage and Transport Conditions: Temperature range: -20 °C to +60 °C / -4 °F to +140 °F		Single Use Only
	Fragile / Handle With Care		UN 2014 – Hydrogen Peroxide
	Do not disassemble		Dangerous Voltage
	Separate collection for electrical and electronic equipment.		Keep Dry
	Keep Out of Direct Sunlight		This Way Up
	Batch Number		Product Number
	Serial Number		Expires (year and month)
	Legal Manufacturer		Date of Manufacture
	Oxidizer – 5.1		Corrosive – 8
	Warning: Hot Surface		Warning: Moving parts, do not touch mechanism
	Cannot be transported by air freight		Wear Gloves

	Environmental Conditions: Operating temperature range for the trophon device: 17 °C to 27 °C / 62.6 °F to 80.6 °F		AcuTrace® RFID Zone
	EU Importer	EC REP	European Authorised Representative
	Conforming to Medical Devices Regulation (EU)2017/745 (MDR)		Conforming to RoHS 3 (EU 2015/863)
	Medical Device		
Touch Screen Symbols			
	Start up from sleep		Cycle Start
	Connected to network		Disconnected from network
	trophon software download		Menu
	No Data Transmission to PACS Server		
   Integrated Probe Positioning Guide			

The following section of the indications for use has been updated relating to the Minimum Operational Cycle Time.

A2.1 Indications for Use

The trophon2 device is intended for the high-level disinfection (HLD) of non-lumened, reusable, transiently invasive and non-invasive medical instruments/devices* e.g. devices that are intended for use for imaging, diagnostic, ablation, coagulation and their accessories.

The trophon2 system consists of a multiple-use device, combined with a single-use disinfectant “trophon NanoNebulant”, delivered from a multi-dose cartridge.

The trophon2 device is suitable for use in general hospital and health care facilities by trained personnel.

The trophon NanoNebulant should be used with the following contact conditions:

Minimum Operational Cycle Time: 4 minutes

Minimum Concentration: 31.5%

Minimum Disinfectant Dose: 1.0 g

Minimum Chamber Temperature: 56 °C

The trophon2 device is NOT intended to reprocess single-use probes or instruments, or to pre-clean medical instruments.

Chemical indicator use is required with every HLD cycle. Only the trophon Chemical Indicator product is the approved chemical indicator for use with trophon2 device.

* The terms “ultrasound probe” and “probe” in the User Manual refer to approved medical instruments/devices.

Part B – SETUP

The following Installation Guide section of the user manual has been added for a new feature related to DICOM.

B2.6 Basic Settings Network

For DICOM user, follow additional instructions as below:

- **DICOM**

To view DICOM network settings, select:

Menu → Settings -> Network -> DICOM

When a trophon2 Plus Software Upgrade device is not connected to a DICOM PACS server, it is indicated at the top right of the idle screens, by the following symbol:

To configure DICOM server list, select:

Menu → Settings -> Network -> DICOM -> Servers

To configure DICOM server settings, select:

Menu → Settings -> Network -> DICOM -> Services

To check subscription and other DICOM settings, select:

Menu → Settings -> Network -> DICOM -> Setting

Part C – CYBERSECURITY

Protect your hospital network and restrict unauthorized access to both the network and any connected medical devices. Monitor network activity and keep antivirus and firewalls updated.

Consult Nanosonics Service if you notice any inconsistencies or strange behaviour from the trophon2 Device.


It is recommended to connect the trophon2 Device to the hospital network via the ethernet, to ensure software updates can be made. Software updates are critical to ensure the trophon2 Device has the latest features and best protection available.

Consult your IT department to set up the network at time of installation. If required, contact your customer service representative to request a copy of the Manufacture Disclosure Statement for Medical Device Security (MDS2).

Trophon2 Cybersecurity Architecture

The Trophon2 Cyber Architecture (Figure 1) gives a pictorial view of the connections made by and to the Trophon2 device in a network environment.

Figure 1: Cybersecurity Architecture

Trophon2 SBOM

If required, contact your customer service representative to request SBOM information.

Trophon2 Software Updates

Trophon2 Device Software Update Methods:

1. **Authorized Service Technician Update:**
 - A certified technician performs the software update on-site.
2. **Remote Software Update:**
 - When connected to the Nanosonics cloud service, the Trophon2 device can receive updates remotely.
 - Users are notified via the device GUI when a new update is available.
 - Users can choose to proceed with or disregard the update.

Events and logging

- **Key Device Behaviors:** Logs include setting changes, operations running an HLD, and etc.
- **Abnormal Events:** Recorded in device logs, accessible only by a Nanosonics representative using specialized software on an external device.
- **Cloud Connection Activity:** Logged in the device internal logs, available via authenticated connection.

Trophon2 Device protective features

The trophon2 Device implements strong cryptographic protection on the network communications. Network communications are limited to only the essential web services for the Service and Maintenance (S&M) software. Those S&M services require approved security certificates to allow access. All other network communications are ignored.

Backup and Restore features / restoring authenticated user

Your trophon2 Device configuration is verified along with the device software to create an intrinsic set of parameters used to operate the device. Since the device configuration is part of the software, there is no need to backup, restore or change this configuration.

Log Files

Any abnormal events will be recorded in the device logs. These are not user accessible but may be accessed by a Nanosonics representative using appropriate software on an external device.

Trophon2 Device end of life

The trophon2 Device software is prepared by Nanosonics and does not depend upon external components for operation. As such, Nanosonics can ensure continued support for upgrades of the device software into the future.

Secure decommission of the trophon2 Device

There is no patient information stored on the trophon2 Device. Decommissioning of the device is acceptable without erasure of data.

Trophon2 Device Software Licenses

The trophon2 device software uses statically linked components that are under MIT, GPL-2.0, LGPL 2.1, BSD 3-clause, EPL 1.0, GPL-3.0 Apache-2.0 licenses.

- MIT license [1]: Yocto Project-Dunfell, libudev
- GPL 2.0 license [2]: Linux Kernel, U-Boot
- Apache-2.0 license [3]: Libraries used in Java and Qt backend, eg; spring, grpc, netty, httpclient, dbus, quartz, log4j, derby etc
- BSD 3-clause [4]: protobuf-java, libupb
- LGPL 2.1 [5]: hibernate, libz
- EPL1.0 [6]: junit
- GPL 3.0 [7]: libgpr

Please see below for links to the referenced license notes:

- [1] <https://opensource.org/license/mit>
- [2] <https://www.gnu.org/licenses/old-licenses/gpl-2.0.txt>
- [3] <https://www.apache.org/licenses/LICENSE-2.0.txt>
- [4] <https://opensource.org/license/bsd-3-clause>
- [5] <https://www.gnu.org/licenses/old-licenses/lgpl-2.1.en.html>
- [6] <https://www.eclipse.org/legal/epl/epl-v10.html>
- [7] <https://www.gnu.org/licenses/gpl-3.0.en.html>

GLOSSARY

PACS

Picture Archiving and Communication System

DICOM

Digital Imaging and Communications in Medicine

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**Manufactured by:**

Nanosonics Limited
7-11 Talavera Road,
Macquarie Park NSW 2113,
Australia
+61 2 8063 1600
customerservice@nanosonics.com.au
www.nanosonics.com.au

**Authorised European Representative/ Importer**

Nanosonics Europe GmbH
Poppenbütteler Bogen 66
22399 Hamburg
Germany
+49 40 46856885
customerservice@nanosonics.eu
www.nanosonics.eu

**UK Responsible Person:**

Nanosonics UK Limited
Ground Floor at The Forum
Unit C1 & C2, Hercules Business Park,
Bird Hall Lane, Stockport, SK3 0UX, UK
+44 (0) 161 686 3030
customerservice@nanosonics.co.uk
www.nanosonics.co.uk

Conforming to WEEE Directive 2012/19/EU under Article 7 Recovery

www.nanosonics.com