trophon[®]2

User Manual Addendum

for trophon[®]2 Plus

Software Upgrade

ENGLISH.....	3
FRANÇAIS	23
ESPAÑOL.....	43

trophon[®]2

User Manual Addendum

for trophon[®]2 Plus

Software Upgrade

Read this addendum together with trophon®2 user manual before operating device to determine the correct procedures.

The Table of Contents in this addendum lists updates to the trophon®2 user manual.

For further information, contact your customer service representative or visit the Nanosonics website.

© 2025 Nanosonics Limited. All rights reserved.

www.nanosonics.com


Table of Contents

Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS	6
A1.1 Labels and Symbols.....	6
A2.1 Indications for Use	8
Part B – SETUP	9
B2.3 Initial Setup	9
B2.6 Basic Settings Network	9
Part E – MAINTENANCE AND ROUTINE CARE.....	10
SECTION E4: Cleaning and Disinfecting.....	10
Part H – CYBERSECURITY	11
Part I – EMC Compliance	14
APPENDIX 1: trophon2 device Technical Specifications.....	21
GLOSSARY	22

Part A – WARNINGS, INTRODUCTION AND INSTRUCTIONS

The following section of the labels and symbols table has been updated to add a new feature relating to DICOM in Touch Screen Symbols.

A1.1 Labels and Symbols

	Caution		Warning
	Consult Instructions For Use		Corrosive
	Environmental Conditions: trophon Storage and Transport Conditions: Temperature range: -20 °C to +60 °C / -4 °F to +140 °F		Single Use Only
	Fragile / Handle With Care		UN 2014 – Hydrogen Peroxide
	Do not disassemble		Dangerous Voltage
	Separate collection for electrical and electronic equipment.		Keep Dry
	Keep Out of Direct Sunlight		This Way Up
	Batch Number		Product Number
	Serial Number		Expires (year and month)
	Legal Manufacturer		Date of Manufacture
	Oxidizer – 5.1		Corrosive – 8
	Warning: Hot Surface		Warning: Moving parts, do not touch mechanism
	Cannot be transported by air freight		Wear Gloves

	Environmental Conditions: Operating temperature range for the trophon device: 17 °C to 27 °C / 62.6 °F to 80.6 °F		AcuTrace® RFID Zone
Touch Screen Symbols			
	Start up from sleep		Cycle Start
	Connected to network		Disconnected from network
	trophon software download		Menu
	No Data Transmission to PACS Server		
1 2 3 Integrated Probe Positioning Guide			

The following section of the indications for use has been updated relating to the Minimum Operational Cycle Time.

A2.1 Indications for Use

The trophon2 is designed to provide High-Level Disinfection (HLD) of validated ultrasound probes. High-Level Disinfection is achieved by surface exposure to a controlled dose of hydrogen peroxide mist delivered to a disinfection chamber contain the ultrasound probe.

The trophon2 system consists of a multiple use instrument combined with a single use disinfectant “trophon Sonex-HL”, delivered from a multi-dose cartridge.

The trophon2 is suitable for use in general hospital and health care facilities by trained personnel.

The trophon Sonex-HL should be used with the following contact conditions:

Minimum Operational Cycle Time: 4 minutes

Minimum Concentration: 31.5%

Minimum Disinfectant Dose: 1.0 g

Minimum Chamber Temperature: 56 °C

Part B – SETUP

The following Installation Guide section of the user manual has been added for a new feature related to DICOM.

B2.3 Initial Setup

To access the trophon2 device software version and build number:

Select **Menu -> Information -> Device**

B2.6 Basic Settings Network

For DICOM user, follow additional instructions as below:

- **DICOM**

To view DICOM network settings, select:

Menu → Settings -> Network -> DICOM

When a trophon2 Plus Software Upgrade device is not connected to a DICOM PACS server, it is indicated at the top right of the idle screens, by the following symbol:

To configure DICOM server list, select:

Menu → Settings -> Network -> DICOM -> Servers

To configure DICOM server settings, select:

Menu → Settings -> Network -> DICOM -> Services

To check subscription and other DICOM settings, select:

Menu → Settings -> Network -> DICOM -> Setting

To clear DICOM cache records and restore default DICOM configuration, select:

Menu → Settings -> Network -> DICOM -> Maintenance

To enable or disable DICOM traceability to the Nanosonics cloud service, select:

Menu → Settings -> Network -> DICOM -> Traceability Data

Part E – MAINTENANCE AND ROUTINE CARE

SECTION E4: Cleaning and Disinfecting

- Do NOT submerge the trophon2 device, or pour liquids over it.
- Keep the trophon2 device level and upright at all times.
- Keep the power socket completely dry (see Figure 4).

Ensure appropriate PPE, including clean gloves, is worn when cleaning and disinfecting the trophon2 device.

For cleaning, Nanosonics recommends that the surfaces of the trophon2 device be thoroughly cleaned in accordance with organisational policies and following the detection of visible contamination. When cool, thoroughly wipe the chamber including the silicone seals, IPP and cable clamp and outside surfaces of the trophon2 device including the screen, handle and underneath the cable tray (refer to section B1.2 Cable Tray for tray removal) with a cloth or wipe moistened with a mild general purpose detergent, or soap and water solution, such as dish soap or general cleaners, until all surfaces are visibly clean. The cleaning process should be repeated if the trophon2 device is not visually clean. Refer to the product's IFU for general instructions and directions for use.

For disinfecting, Nanosonics recommends that the surfaces of the trophon2 device be disinfected in accordance with organisational policies. Wipe all accessible surfaces of the trophon2 device with an Isopropanol or Quaternary Ammonium (Quat) wipe. Refer to the wipe product's IFU for general instructions and directions for use. Ensure surfaces remain wet for the full contact time for disinfection according to the wipe manufacturer's IFU and use additional wipes if required.

Part H – CYBERSECURITY

The trophon2 device has built-in cybersecurity features designed to help protect both the device and your hospital network. To maintain the highest level of protection, ensure the device is updated to the latest software version.

It is recommended to connect the trophon2 Device to the hospital network via the ethernet, WiFi or Cellular to ensure timely software updates that provide the latest features and security enhancements.

Consult your IT department to to network your device at time of installation. If required, you can request a copy of the Manufacture Disclosure Statement for Medical Device Security (MDS2) from your customer service team at customerservice@nanosonics.us.

trophon2 Cybersecurity Architecture

The trophon2 Cyber Architecture (Figure 1) gives a pictorial view of the connections made by and to the trophon2 device in a network environment.


Figure 1: trophon2 Cybersecurity Architecture

trophon2 Software Updates

trophon2 device software update methods:

- Software updates can be applied either by a certified service technician on-site, via self-serve tool, or remotely when the device is connected to the Nanosonics cloud service.
- When a new update is available, the device notifies the user via the touchscreen. Users can choose whether to proceed with the update.

System Monitoring and Log Management

The trophon2 device automatically records operational activity and any abnormal events in internal logs. These logs are accessible only to Nanosonics representatives using specialized tools.

Consult Nanosonics Service if you notice any inconsistencies or strange behaviour from the trophon2 Device

trophon2 Device protective features

Communications between the device and approved service applications are encrypted and restricted to essential functions. Only connections with valid security certificates are accepted. All other network traffic is blocked.

Backup and Restore features / restoring authenticated user

Your trophon2 Device configuration is verified along with the device software to create an intrinsic set of parameters used to operate the device. Since the device configuration is part of the software, there is no need to backup, restore or change this configuration.

Log Files

Any abnormal events will be recorded in the device logs. These are not user accessible but may be accessed by a Nanosonics representative using appropriate software on an external device.

trophon2 Device end of life

The trophon2 Device software is prepared by Nanosonics and does not depend upon external components for operation. As such, Nanosonics can ensure continued support for upgrades of the device software into the future.

Secure decommission of the trophon2 Device

There is no patient information stored on the trophon2 Device. Decommissioning of the device is acceptable without erasure of data.

Additional Information

A copy of Software Bill of Materials (SBOM), Manufacture Disclosure Statement for Medical Device Security (MDS2), and other software license details can be provided on request by contacting Nanosonics Customer Service at customerservice@nanosonics.us

trophon2 Device Software Licenses

The trophon2 device software uses statically linked components that are under MIT, GPL-2.0, LGPL 2.1, BSD 3-clause, EPL 1.0, GPL-3.0 Apache-2.0 licenses.

- MIT license [1]: Yocto Project-Dunfell, libudev
- GPL 2.0 license [2]: Linux Kernel, U-Boot

- Apache-2.0 license [3]: Libraries used in Java and Qt backend, eg; spring, grpc, netty, httpclient, dbus, quartz, log4j, derby etc
- BSD 3-clause [4]: protobuf-java, libupb
- LGPL 2.1 [5]: hibernate, libz
- EPL 1.0 [6]: junit
- GPL 3.0 [7]: libgpr

Please see below for links to the referenced license notes:

- [1] <https://opensource.org/license/mit>
- [2] <https://www.gnu.org/licenses/old-licenses/gpl-2.0.txt>
- [3] <https://www.apache.org/licenses/LICENSE-2.0.txt>
- [4] <https://opensource.org/license/bsd-3-clause>
- [5] <https://www.gnu.org/licenses/old-licenses/lgpl-2.1.en.html>
- [6] <https://www.eclipse.org/legal/epl/epl-v10.html>
- [7] <https://www.gnu.org/licenses/gpl-3.0.en.html>

Part I – EMC Compliance

FCC ID: 2AM5R-TROPHON2

IC: 22999-TROPHON2

FCC RF Exposure Statement

This equipment complies with the FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 cm (8 inches) between the radiator and all persons during normal operation.

Notice

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Warning

Any changes or modifications not expressly approved by Nanosonics could void the user's authority to operate this equipment.

The trophon2 device includes a HF RFID reader that operates at 13.56MHz.

Complies with:

IEC 60601-1-2:2020

IEC TS 60601-4-2:2024

IEC 61326-1:2020

CISPR 11:2024

FCC Part 15 Subpart C Section 15.225

47 CFR Part 15 Subpart B

47 CFR Part 2.1091

This device is intended to be used in professional healthcare facility environments only. It must not be used on aircraft. It is intended to automatically perform High Level Disinfection (HLD) of the validated ultrasound probes. If you suspect that trophon2 is being affected by other devices located close by, try moving the other device further away. If the issue persists then contact service support.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection from radio-frequency communication services. The user may need to take mitigation measures, such as relocating or re-orienting the equipment.

Warning

Use of the trophon2 device adjacent to, or stacked with, other equipment should be avoided as it may result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of the trophon2 device may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

The trophon2 device is only to be used with the supplied mains cord and a Cat 5e or better Ethernet cable.

Warning

Portable RF communications equipment (including peripherals such as antenna cables and external antennae) should be used no closer than 30 cm (12 inches) to any part of the trophon2 device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic Emission Compliance Information and Guidance

Guidance and Manufacturer Declaration – Electromagnetic Emissions		
The trophon2 device and related components are intended for use in the following prescribed electromagnetic environments only. Customers and users of this product are advised to check that it is being used in such environments.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
Radiated Emissions CISPR 11	Class A (Group 1)	The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals only. It is not suitable for use in domestic environments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Conducted Emissions CISPR 11	Class A (Group 1)	
Harmonic Emissions IEC 61000-3-2	Not Applicable	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Not Applicable	

The trophon2 device includes a HF RFID reader that operates at 13.56MHz using On-off Keying (OOK). The maximum Effective Radiated Power (ERP) is 1.1µW.

Electromagnetic Immunity Compliance Information and Guidance

The **Essential Performance** of the trophon2 is defined as being free from indicating a successfully completed HLD cycle when such a cycle was not completed.

Guidance and Manufacturer Declaration – Electromagnetic Immunity			
The trophon2 device and related components are intended for use in the following prescribed electromagnetic environments only. Customers and users of this product are advised to check that it is being used in such environments.			
Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV air ± 4 kV air ± 8 kV air ± 15 kV air	± 8 kV contact ± 2 kV air ± 4 kV air ± 8 kV air ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. In the event of the Ethernet connection being lost please power cycle the trophon2
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. In the event of the printer going offline please power cycle the trophon2
Surge IEC 61000-4-5	0.5 and 1.0 kV (Line to Line) 0.5, 1.0 and 2.0 kV (Line to Earth) 2.0 kV for Signal Input/Output (Line to Earth)	0.5 and 1.0 kV (Line to Line) 0.5, 1.0 and 2.0 kV (Line to Earth) 2.0 kV for Signal Input/Output (Line to Earth)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/ Dropout IEC 61000-4-11	0% UT for 0.5 cycles (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) @ 50/60 Hz 0% UT for 1 cycle (0°) @ 50 /60 Hz 70% UT (30% dip in UT) for 25/30 cycles (0°) @ 50 /60 Hz	0% UT for 0.5 cycles (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) @ 50/60 Hz 0% UT for 1 cycle (0°) @ 50 /60 Hz 70% UT (30% dip in UT) for 25/30 cycles (0°) @ 50 /60 Hz	Mains power quality should be that of a typical commercial or hospital environment.

	0% UT for 250/300 cycles (0°) @ 50 /60 Hz	0% UT for 250/300 cycles (0°) @ 50 /60 Hz	
Radiated RF IEC 61000-4-3	80 MHz to 2.7 GHz 3 V/m 380 – 290 MHz 27 V/m 430 – 470 MHz 28 V/m 704 – 787 MHz 9 V/m 800 – 960 MHz 28 V/m 1700 – 1990 MHz 28 V/m 2400 – 2570 MHz 28 V/m 5100 – 5800 MHz 9 V/m	80 MHz to 2.7 GHz 3 V/m 380 – 290 MHz 27 V/m 430 – 470 MHz * 28 V/m 704 – 787 MHz * 9 V/m 800 – 960 MHz 28 V/m 1700 – 1990 MHz * 28 V/m 2400 – 2570 MHz 28 V/m 5100 – 5800 MHz 9 V/m	Portable RF communications equipment (including peripherals such as antenna cables and external antennae) should be used no closer than 30 cm (12 inches) to any part of the trophon2 device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:  *At these frequencies and field strengths above 3 V/m, the device may display an error message and indicate the disinfection cycle has failed. If this issue reoccurs, move any RF transmitters further away from the device.
Conducted RF IEC 61000-4-6	0.15 to 80 MHz, 3V rms 6 V in ISM bands between 0.15 MHz and 80 MHz	0.15 to 80 MHz, 3V rms 6 V in ISM bands between 0.15 MHz and 80 MHz	
Proximity magnetic fields IEC 61000-4-39	134.2 kHz, Pulse modulation 2.1 kHz 65 A/m 13.56 MHz, Pulse modulation 50 kHz 7.5 A/m	134.2 kHz, Pulse modulation 2.1 kHz 65 A/m 13.56 MHz, Pulse modulation 50 kHz 7.5 A/m	Suitable for use in a typical commercial or hospital environment. External RFID readers should not be placed next to the device. An error message may be displayed and indicate the disinfection cycle has failed. If this issue reoccurs, move any RFID readers further away from the device.

The **Intended Use** of the trophon2 is to provide high-level disinfection (HLD) of validated ultrasound probes. Other functions include the scanning of RFID tags for traceability and record keeping and the option of connecting USB peripherals. There is also the option of connecting trophon2 to a hospital network using ethernet.

Guidance and Manufacturer Declaration – Electromagnetic Immunity For Intended Use			
The trophon2 device and related components are intended for use in the following prescribed electromagnetic environments only. Customers and users of this product are advised to check that it is being used in such environments.			
Immunity Test	IEC TS 60601-4-2 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 4 kV contact ± 2 kV air ± 4 kV air ± 8 kV air	± 4 kV contact ± 2 kV air ± 4 kV air ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	0.5 and 1.0 kV (Line to Line) 0.5, 1.0 and 2.0 kV (Line to Earth) 2.0 kV for Signal Input/Output (Line to Earth)	0.5 and 1.0 kV (Line to Line) 0.5, 1.0 and 2.0 kV (Line to Earth) 2.0 kV for Signal Input/Output (Line to Earth)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/ Dropout IEC 61000-4-11	0% UT for 0.5 cycles (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) @ 50/60 Hz 0% UT for 1 cycle (0°) @ 50 /60 Hz 70% UT (30% dip in UT) for 25/30 cycles (0°) @ 50 /60 Hz 0% UT for 250/300 cycles (0°) @ 50 /60 Hz	0% UT for 0.5 cycles (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) @ 50/60 Hz 0% UT for 1 cycle (0°) @ 50 /60 Hz 70% UT (30% dip in UT) for 25/30 cycles (0°) @ 50 /60 Hz 0% UT for 250/300 cycles (0°) @ 50 /60 Hz	Mains power quality should be that of a typical commercial or hospital environment.

Radiated RF IEC 61000-4-3	80 MHz to 6 GHz 3 V/m	80 MHz to 6 GHz 3 V/m	Portable RF communications equipment (including peripherals such as antenna cables and external antennae) should be used no closer than 30 cm (12 inches) to any part of the trophon2 device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	0.15 to 80 MHz, 3V rms	0.15 to 80 MHz, 3V rms	

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the trophon2 device.

The trophon2 device and related components are intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The Customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the trophon2 device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (in metres)		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

APPENDIX 1: trophon2 Device Technical Specifications

Corrected typographical error in electromagnetic compliance section for product class.

Electrical Specification	Rated mains input voltage: 120 V AC Rated mains input current: 6 A, 50/60 Hz Mains Inlet: IEC type C13 Equipment must be connected to an earthed outlet using the power cable supplied with the trophon2 device.
Data port	Ethernet connector RJ45 USB Port: Type A
Environmental Specification	Operating temperature range: 17 °C to 27 °C / 62.6 °F to 80.6 °F
Storage and Transport Conditions	Temperature range: -20 °C to +60 °C / -4 °F to +140 °F
Physical Characteristics	Weight of trophon2 device: Unpacked: 22 kg (48.5 lbs) Packed: 25 kg (55 lbs) Dimensions of trophon2 device: 535 mm high × 360 mm wide × 317 mm depth (21 in high × 14.2 in wide × 12.5 in depth)
Electromagnetic Compliance	The trophon2 device has been tested and found to comply with the limits for emission (electromagnetic Interference) pursuant to EN61326-1:2013 (CISPR 11 Group 1 Class A limits)

GLOSSARY

DICOM

Digital Imaging and Communications in Medicine

EMC

Electromagnetic Compatibility

PACS

Picture Archiving and Communication System

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Nanosonics Inc.